AMENDMENT

IN THE CLAIMS:

Please amend the claims as follows:

- (Currently Amended) A method for detecting pathogenic mycobacteria in clinical specimens, said method comprising the steps of:
 - (a) clarifying the clinical specimens from eontainment <u>contaminant</u> by conventional methods,
 - (b) treating the processed clinical specimens obtained in step (a) with a modified lysis buffer to inactivate live pathogenic mycobacteria to make the process safe for the user.
 - (c) extracting genomic DNA from the processed clinical specimen obtained from step (b) using a modified method to increase the yield and quality of DNA,
 - (d) designing sequence of SEQ ID NO: 4 from the DNA obtained in step (c) for specific detection of pathogenic mycobacteria, said designed sequence comprising of selected intergenic region of SEQ ID. NO: 3, a flanking region containing a portion of the gene mmaA1 of SEQ ID NO: 1 and a portion of gene mmaA2 of SEQ ID NO: 2 of the DNA obtained in step (c).
 - (e) designing and synthesizing a set of specific oligonucleotide primers of SEQ ID NO: 5, which is the forward primer and SEQ ID NO: 6, which is the reverse primer for Polymerase Chain Reaction (PCR) amplification of SEQ ID NO: 4,
 - (f) developing a PCR amplification process for specific amplification of SEQ ID NO: 4 of (d), said process comprising using the specific oligonucleotide primers designed and synthesized in step (e) for detecting presence of pathogenic mycobacteria in the clinical specimens, and
 - (g) analyzing the amplified PCR product by restriction fragment length polymorphism (RFLP) analysis for differentiation of the species of the pathogenic mycobacterium for a quick assessment of HIV co-infection.

2. (Previously Presented) A method as claimed in claim 1, wherein the designed SEQ ID NO: 4 has a sequence as follows:

5'TGGATCCGTTGACCATGAGGTGTAATGCCTTTCCGGACCCTAGGTGGCCTTT
CGGTGCTTGCACGGAACGCACCGATGCTTCCCCCCTCCCCGCATGCTCGAGGC
ATGCTATCCGATACAGGGCCGCCGACTAAACCGCGATCGAATTTGCCCAGG
TCAGGGAACGGATATGAGCGGACGAGCTACTTGGTCATGGTGAACTGGGCGA
CGTTGATTAGGCCTCTGCGGAAGCGCTCCGCGCATCCGGTCAGATAGTGCAT
GAAGTTGTTGTAGACCTCTTCGGACTGTACGGCGATGGCGCGTTCGCGGGCA
GCCTGTAGGTTGGCGGCCCATGCATCGAGAGTCCGTGCGTAGTGGGAATTC
3'.

- 3. (Previously Presented) A method as claimed in claim 1, wherein the clinical specimen is selected from the group consisting of sputum, gastric lavage, cerebrospinal fluid, blood, tissue biopsies, bone marrow aspirates and other body fluids or tissues.
- 4. (Previously Presented) A method as claimed in claim 1, wherein clarification of the specimens in step (a) from the contaminants is carried out by adding to said specimens a digestion decontamination mix containing mild alkali, NaOH, tri sodium citrate and a mucolytic agent and guanidinium isothiocyanate in the range of about 0.4-2.5 M followed by concentrating the specimens by centrifugation.
- 5. (Previously Presented) A method as claimed in claim 4, wherein the digestion decontamination mix containing mild alkali, NaOH, tri sodium citrate and a mucolytic agent and guanidinium isothiocyanate in the range of about 0.5-2.0 M.
- 6. (Previously Presented) A method as claimed in claim 1, wherein the DNA in step (c) is extracted from the treated clinical specimen using a modified lysis buffer by inclusion of ingredients comprising guanidinium isothiocyanate in a range of about 0.5-8 M, Tris.Cl pH 7.6 in

a range of about 20-100 mM, N lauryl Sarcosyl in a range of about 0.5-2% by weight of the buffer, EDTA in a range of about 0.1-20 mM, β -Mercaptoethanol in a range of about 1-25 mM and NaCl is present in an amount of about 0.2M; and purifying the DNA to improve yield by thorough precipitation by organic solvents.

- 7. (Previously Presented) A method as claimed in claim 6, wherein guanidinium isothiocyanate is about 4M, Tris-HCl pH 7.6 is about 50 mM, N lauryl Sarcosyl is 1% by weight of the buffer, EDTA 1 mM, β-Mercaptoethanol is about 10 mM and NaCl is about 0.2M.
- 8. (Canceled).
- (Previously Presented) A method as claimed in claim 1, wherein the genomic DNA yield is increased 25 to 50%.
- 10. (Canceled).
- 11. (Previously Presented) A method as claimed in claim 1, wherein high yielding amplification of DNA in step (f) is achieved by the modified Touch Down PCR cycling conditions, said conditions comprising steps of providing an initial high annealing temperature in the range of 62-72°C followed by lowering of temperature in the range of 0.1-1°C per PCR cycle for the first 10-25 cycles, then subsequently carrying out 30 PCR cycles at an optimum annealing temperature of 56-62°C.
- 12. (Previously Presented) A method as claimed in claim 1, wherein high yielding amplification of DNA is achieved by modified Touch Down PCR cycling conditions, said conditions comprising steps of providing an initial high annealing temperature of 70°C followed by lowering of temperature of 0.8°C per PCR cycle for about first 14 cycles to about 58°C for another 25 PCR cycles.

PATENT USSN: 10/725,994

Atty Dkt: 041144.010

13. (Previously Presented) A method as claimed in claim 1, wherein the oligonucleotide primers capable of amplification of intergenic region of SEQ ID NO: 4 for detection of pathogenic Mycobacteria in clinical specimens are selected from group consisting of:

- a. 5' TGGATCCGTTGACCATGAGGTGTAATG 3' (SEQ ID NO: 5), which is the forward primer, and
- b. 5' GGAATTCCACTACGCACGGACTCTC 3' (SEQ ID NO: 6), which is the reverse primer.
- 14. (Previously Presented) A method as claimed in claim 1, wherein the length of oligomeric primers is between 5 and 100 bases.
- 15. (Previously Presented) A method as claimed in claim 1, wherein the modified lysis buffer provides a cleaner preparation of the DNA.
- 16. (Previously Presented) A method as claimed in claim 1, wherein treatment with the modified lysis buffer containing 4M guanidinium isothiocyanate inactivates the live mycobacteria to make the procedure safer for the operator.
- 17. (Currently Amended) A diagnostic kit for the detection of pathogenic mycobacteria in clinical specimens, comprising primers selected from the group-consisting of the following primers:
 - a. 5' TGGATCCGTTGACCATGAGGTGTAATG 3' (SEQ ID NO: 5), which is the forward primer, and
 - b. 5' GGAATTCCACTACGCACGGACTCTC 3' (SEQ ID NO: 6), which is the reverse primer.
- 18. (Previously Presented) A method as claimed in claim 1, wherein the contaminant clarified in step (a) comprises mucus and/or live organisms other than mycobacteria.

(Currently Amended) A set of primers of SEQ ID NOs: 5 and 6 eemprising consisting of:
 5' TGGATCCGTTGACCATGAGGTGTAATG 3' (SEQ ID NO: 5), which is forward

primer; and

 5° GGAATTCCACTACGCACGGACTCTC 3° (SEQ ID NO: 6), which is the reverse primer.